

# Mash Cleaning

Ryan Duve


November 3, 2014

## Abstract

Only workers highly familiar with the technical aspects of Hifrost are authorized to perform this procedure.

We suspect the dilution cooldown (10/9/2014) resulted in 10 L air getting stuck in the recovered mash in T1. This procedure serves as a technical document to capture the air in a LN<sub>2</sub> trap and measure the captured air in the tank T3.


## 1 Prep




Step 1: Red and blue volumes need to be evacuated ahead of time, blue connections are new and must be leak checked via HVLD (“hand valve - leak detector”).

- ☐ Regenerate LN<sub>2</sub> trap.
- ☐ Configure the system as shown in Step 1. Purge all lines along the red/blue paths with P7 to clean the system.
- ☐ Leak check the blue volumes.
- ☐ Close all valves; kill P7; fill trap with LN<sub>2</sub>.


## 2 Cleaning




Step 2: Open valves in this order: HV1, HV19, HV20, HV3, HV4, HV5, HV14, HV13, HV12, HV9, HVBB. Wait until Big Blue equilibrates with T3. The green volumes should all be at the same pressure.



Step 3: Close HV19, HV20; then open HV17 and turn on P10. Slowly open HV19 and let M5 drop.



Step 4: When M5 is around 0 mbar, close HV19, then kill P10. Close HV17. Open HV20 and HV19. Wait until M5 and MBB equilibrate.



Step 5: Close HV19 and HV20. Open HV18 and then start P10. Slowly open HV20; wait until MBB approaches 0. Close HV20, kill P10, then close HV18.

